


Inion CPSTM Baby System

Biodegradable Fixation System

Inion CPS™ Baby System

All Inion CPS™ implants are composed of Inion® biodegradable polymer blends. They have excellent handling properties, strength and degradation characteristics that support a more natural healing process. Inion CPS™ Baby System is designed especially for the use in paediatric craniofacial trauma and reconstruction procedures.

Materials




The implants are composed of L-Lactide, Polyglycolide and TMC (Trimethylene Carbonate).

The proportion of each polymer varies according to the intended application of the specific implant, so that the strength, malleability and degradation profile best suits the clinical requirements. In the case of Inion CPSTM Baby products, degradation is designed to occur faster to prevent growth restriction of the skull bones.

Degradation

Inion CPSTM Baby implants retain minimum of 70% of their initial strength 6 weeks after implantation. Bioresorption takes place within two to three years.



In vitro in phosphate buffer solution at 37°C

Advantages

Inion CPS™ implants have more than 14 years history of safe clinical use, and more than 60,000 operations have been conducted with these products.



Cranial vault remodelling for craniosynostosis










Cranial reconstruction after trauma









Cranial reconstruction after tumor resection









Congenital and acquired cranial deformity









Autologous cranioplasty









Courtesy of Dr. Dylan Murray, Children's University Hospital, Temple Street, Dublin, Ireland

Ordering Information


Plates

Art. No. Description PLT-1000 4 hole plate

PLT-1001 20 hole plate (reinforced centre hole)

PLT-1029 Mesh plate, 14 x 14 holes

PLT-1059 64 hole plate



Screws

Art. No. Description

SCR-1220 1.5 \times 4 mm Screw ring 5 pcs (+ one 2.0 \times 5 mm emerg.) SCR-1221 1.5 \times 6 mm Screw ring 5 pcs (+ one 2.0 \times 7 mm emerg.)

SCR-1210 2.0×5 mm Screw ring 5 pcs SCR-1211 2.0×7 mm Screw ring 5 pcs



Instruments

Art. No. Description

1.5 mm

INS-9116 1.2 mm drill bit short (J-latch)

INS-9002 1.2 mm drill bit with 5 mm stop (J-latch)

INS-9027 1.5 mm bone tap (manual)

INS-9261 1.5 mm self-drilling bone tap (pentagonal) INS-9263 1.5 mm self-drilling bone tap (manual)

2.0 mm

INS-9001 1.6 mm drill bit short (J-latch)
INS-9030 2.0 mm bone tap (manual)

INS-9262 2.0 mm self-drilling bone tap (pentagonal) INS-9264 2.0 mm self-drilling bone tap (manual)



INS-9029 Universal screw driver blade INS-9007 Universal screw driver handle

INS-9024 Plate bending pliers INS-9047 Counter sink 1.5/2.0 mm

ACC-9801 Inion Thermo™ (water bath 230 V)
ACC-9804 Inion Thermo™ (water bath 110 V)
ACC-9802 Inion Thermo™ drape (5/package)
ACC-9813 Inion universal instrument tray
ACC-9818 Inion compact instrument tray
PDSD-5000-1 Pro-Driver, Pro-Dex surgical driver
PDBP-001-2 Pro-Dex Single Pouch Battery pack x 1









INION OY

Lääkärinkatu 2 FI-33520 Tampere, FINLAND tel: +358-10-830 6600

fax: +358-10-830 6601 email: info@inion.com internet: www.inion.com

INION INC

2800 Glades Circle Suite 138, Weston FL 33327

Toll-free tel: 866-INION-US tel: 954-659-9224-fax: 954-659-7997